

Micron advanced imaging course


Ian Dobbie
micron@bioch.ox.ac.uk

Why microscopy?

What can you do with a microscope?

100 years ago:

Magnify small things to visualise more details




Now:

Image the distribution of specific molecules
inside cells


Follow changes in the distribution of molecules
or the morphology of cells over time

Determine how close molecules are to each
other within cells



Useful size range for light microscopy

We use the light microscope to image structures and substructures within this range: from about **300 μm** down to about **0.2 μm (200 nm)**



What is important?

Whats Important in Microscopy ?

Contrast :- What is the difference between what you want to see and everything else?

Resolution :- How small things can you see?


Nothing else

A single point image in a microscope.

A single point imaged on a microscope generates a diffraction pattern

Called a Point Spread Function (PSF)

Central region is the Airy disk, outer bits are Airy rings




Contrast enhanced point spread function

**Why so many different
microscopes?**

FLUORESCENCE MICROSCOPES

Deltavision Elite



New Biochemistry Room 00-031

Widefield fluorescence system

[FIND OUT MORE](#)

Deltavision Core



New Biochemistry Room 00-031

Widefield fluorescence system

[FIND OUT MORE](#)

Spinning Disc (Nasmyth)



New Biochemistry Room 00-031

An Ultraview spinning-disc system equipped with a temperature and CO₂ controlled chamber

[FIND OUT MORE](#)

Spinning Disc (Barr)




New Biochemistry Room 00-031

A fast-acquisition spinning disc system suitable for live or fixed imaging.

[FIND OUT MORE](#)

Zeiss Z1 Lightsheet



New Biochemistry Room 00-060

This is a commercial lightsheet microscope designed primarily for time-lapse imaging of thick specimens.

[FIND OUT MORE](#)

Olympus FV1000



Dunn School of Pathology Room 214.00.25

Olympus FV1000 Laser Scanning Microscope with Becker and Hickel FLIM system.

[FIND OUT MORE](#)

**What is YOUR
application?**

Your application

- The Sample
- Magnification
 - Fixed/live
- Resolution
 - Imaging speed
- Objective
 -
- Fluorescent label